

Amtryke Therapeutic Tricycle Request Form

(Must be filled out completely by adult rider or parent/guardian)

Recipient's Name:	Age:	Date of Request:
Mailing Address:		Phone #:
City/State/Zip:		ail:
Diagnosis:		
If Re	ecipient is Under Age 18	
Parent/Guardian Name:		
If different from above		
Mailing Address:	Pho	ne #:
City/State/Zip:	Email:	
Secondary Contact Name:	Phone #:	
Treating Therapist's Name:		
Phone #:	Email:	
Will you need financial assistance to obtain the left yes, how much can you pay? Note: Amtryke therapeutic tricycles are distributed base tricycles are at the discretion of the local chapter or pare Tell us about the recipient*:	d on available funds and need. Individua	I placements of Amtryke therapeutic
Tell us about the recipient.		
*This information will be made public to help obtain fun	ding. Please don't include information yo	ou don't want shared.
Including a photo of the recipient will help us opereferred but we also accept professionally pring including a photo, you are giving consent for AMBUCS	nted glossy photos. No photocopi	es or folded images.
Signature:		Date:

Individuals will not be considered for placement until all three forms are returned to AMBUCS: this form filled out by the recipient/guardian, as well as the Assessment Form and Tryke Selection Form filled out by the therapist.

Amtryke Therapeutic Tricycle Waiver Form

(Must be filled out completely by adult rider or parent/guardian)

AMBUCS members nationwide are dedicated to creating opportunities for mobility and independence for people with disabilities by providing AmTryke therapeutic tricycles, offering educational scholarships to therapy students and performing various forms of community service.

Purpose:

The AmTryke therapeutic tricycle is designed for people with disabilities. It creates a feeling of freedom, builds self-esteem, strengthens muscles and improves motor coordination and rage of motion—all while making exercise fun.

Steering:

Initially, the rider may have difficulty turning or changing directions. Encourage the rider to go straight ahead, back up and slowly turn around. On many models there are three steering options for the Amtryke. On the front column of the tricycle you will find two holes for the steering pin. The top hole is straight steering, the bottom hole allows a 20 degree turning radius. Leaving the pin out gives the rider free steering.

Safety Cautions

- Fast speeds and sharp turns can cause the Amtryke therapeutic tricycle to tip or turn over.
- Always wear a helmet when riding an Amtryke. Use of other protective gear is highly recommended.
- Adult supervision required if used by younger or developmentally delayed riders.
- Use caution near vehicles, swimming pools and other bodies of water, hills, alleys and sloped driveways.
- Always wear shoes.
- Never allow more than one rider.
- Use of the steering pin is recommended to prevent over-steering and possible tip-overs.

The information contained in this document is not intended nor implied by National AMBUCS $^{\text{TM}}$, Inc. to be professional medical advice by National AMBUCS, Inc. Always seek the advice of your physician, therapist or other qualified healthcare provider prior to staring any treatment or with any question you may have regarding a medical condition. Nothing contained in this document is intended by National AMBUCS, Inc. to be for medical diagnosis or treatment by National AMBUCS, Inc. or on behalf of National AMBUCS, Inc.

In no event shall National AMBUCS, Inc. be liable for any direct, indirect, incidental, consequential, special, exemplary, punitive, or any other monetary or other damages, personal injury or property damages, fees, fines, costs, attorney fees, or liabilities of any kind arising out of or relating in any way to this service or use of the Amtryke® therapeutic tricycle, and/or content or information provided herein.

☐ I agree that I may be photographed by National AMBUCS. I also agree that my photo and name may be used in promo	
efforts for National AMBUCS and Amtryke, LLC. I further grant National AMBUCS the ability to use the photos and name advertising/publicity purposes without additional compensation, except where prohibited by law.	tor
By signing below, I acknowledge that I have read and understood this liability waiver. Recipient's Name (printed):	

Recipient's Name (printed):	
Adult Recipient Signature:	
If Recipient is Under Age 18	
Legal Guardian Name:	
Legal Guardian Signature:	Date:

Individuals will not be considered for placement until all three forms are returned to AMBUCS: this Form filled out by the recipient/guardian, as well as the Assessment Form and Tryke Selection Form filled out by the therapist.

Amtryke Assessment Form (Must be filled out completely by therapist)

Recipient's	s Name:	(*******		completely by	,		
Age:		Weight (lbs.):*		Height (in	ches):*	
Diagnosis(es):*					<i>'</i> –	
		*This private infor	mation is only	used to help ap	propriately fit the ric	der.	
	RIDEF	R'S MEASUREM	ENTS			Α	Center of Shoulder
Arm Mea	asurements (inch	nes)	1	otal Length	Ç	В	Center of Elbow
Left	A to B:	B to C:			A	<u>C</u>	Center of Digit Crease
Right	A to B:	B to C:			-	D E	Center of Hip Center of Knee
Trunk	A to D:					F	Bottom of Foot
Leg Mea	surements (inche	es)	1	otal Length	E	, <u></u>	
Left	D to E:	E to F:					Arm Length & Leg
Right	D to E:	E to F:			∡ F		Length Measurements
							are critical to correct Amtryke Selection
Notes on	Provided Measur	rements (if any):		Helmet S			Sizing Chart is
			Size		ent (head circumferen	ce)	available online:
			XXS		8.5" to 19.5"	,	www.ambucs.org/riders
			XS		20.5" to 22" 22" to 23.6"		/wish-list/sizing-chart/
			S/M L/XL		22 (0 23.6 3.6" to 25.75"		
			L/ //L	2.	5.0 (0.25.75		
Therapist I	Name:			_ Is this the tr	reating therapist?	□ Ye	es 🗆 No
Credential							
Street Add				City:		State:	Zip:
Is this facil		valuation Site?		No □ Not	t Sure		
Therapist	comments conce	erning recipient or g	oals:				
This reque	st/assessment is o	directed to:					
	Local AMBUCS Ch	apter Name:					
		· t (AMBUCS Resource	e Center)				
		•		on this rider wo	uld henefit from an A	Amtryke '	You assume no liability.
Therapist:	-	Tymig that in your pro-	Coolonal Opinio	on this nati wo	ara seriene irom arr	unti y kc.	Date:
			Ship A	mtryke To			
Name/Fa	cility:				Phon	e:	
Street Ad						-	
City:					State):	Zip:

Recipient's Name:	
·	·

Thanks for choosing an Amtryke therapeutic tricycle!

In order to accommodate the widest variety of people and disabilities, Amtryke offers many tryke models and each can be customized in a variety of ways. Following the steps below will help you choose to the perfect tryke for your client from what might seem like a dizzying array of options.

Remember you can always refer to our website, <u>www.amtrykestore.com</u>, or the Amtryke catalogue for more information and product images.

- **Step 1:** Fill out the first page of the Amtryke Assessment Form.
- **Step 2:** Choose the way the tryke will be propelled: **Hand & Foot, Foot,** or **Hand**. Your choice should be based on the rider's ability and therapy goals.

Hand & Foot trykes improve coordination, strength and range of motion. Using all four extremities helps with weakness in any area, even general weakness, and can positively affect tone.

Foot trykes were developed in response to requests from therapists for a traditional tricycle for riders with special needs.

Hand trykes are designed for persons whose lower limbs lack function or those who need special therapy for the upper extremity.

- **Step 3:** Take rider's measurements from the front of this Assessment Form and compare them to our Sizing Chart. (This will narrow the choices considerably)
- **Step 4:** Choose a drive. (If it doesn't mention a choice, then ignore this step.)

Amtrykes come with two drive possibilities: **fixed** drive or **geared** drive. Tryke models have been carefully designed so the drive matches the therapeutic goals of the equipment; therefore all drives are not available on all trykes.

A **fixed** drive, commonly known as a 'fixie,' works on a mechanical level to help individuals make a full pedal rotation. The foot crank is constantly in motion for full therapeutic effect. Coasting is not possible with a positive drive; when limb motion stops, the bike does as well.

Riders of **geared** drive trykes must be able to make a complete pedal rotation on their own. This tryke is suitable for riders who need help with balance and a stable sitting position and who have the cognitive ability to steer successfully and apply the coaster brake or hand brake. The key feature of geared trykes is the ability to coast.

Step 5: Chose any adaptations and/or accessories needed by the rider. Each tryke model can only be customized in the ways noted in its own model section in the Tryke Selection Forms or with the generic accessories listed below. If a customization option is not listed, it is because of design or other practical constraints in offering it on a particular model.

GENERIC ACCESSORIES (not model specific)					
Safety Items	☐ Safety Flag	□ NEW! Medical Cart	☐ Water Bottle w Cage ☐ License Plate		
		(carries 0 ² bottle/ventilator/tubing)			
Leg & Foot Items	☐ Foot Cups (pair):	☐ Knee Adductor Strap:	☐ Pedal Block (1 = ¾")qty		
	□ Small	□ Small □ Large			
	□ Medium	□ Medium			
Hand Items	☐ Variable Range of	□ Wrist Wraps	☐ Wrist Brace Mitt:		
	Motion Kit	(Includes right & left)	☐ Right ☐ Left		
	(only for Hand & Foot Cycles)	□ X-Small □ Large	□ XX-Small □ Medium		
		□ Small	□ X-Small		
		□ Medium	□ Small		

Recipient's Name:			
OOT TRYKES			
□ New! 1410 FOOT	TRYKE – (10" wheels, fixed dr	rive) 50-FC-0100 Arm Leng	th 7-20", Leg Length 14-21"
☐ Option #1 – Standard	Seating System: Blue Bucket	t Seat	
·	Seating System: Snappy Sea		
•	sories: ☐ Notched Seat Bottom	•	□ Head Rest
		m Cushion* □ XL Seat Bacl	
1410 Accessories:	☐ Rear Steering Kit	☐ Pull/Steering	•
	□ Push Bar	☐ Heavy-Duty F	
	YKE (12" wheels, fixed drive) — 5		.3", Leg Length 17-22"
•	l Seating System: Snappy Sea	•	
Snappy Seat Access	sories: □ Notched Seat Bottom		☐ Head Rest
	☐ Notched Wide Botton	m Cushion* □ XL Seat Bac	ck** ☐ Separator Cube
□ Option #2 – Alternate	e Seating System: Blue Bucke	t Seat	
•	e Seating System: Saddle Sea		
Snappy Accessories:	☐ Pull/Steering Bar	,	☐ H-Harness: ☐ 8" ☐ 11.5"
	☐ Push Bar	☐ Heavy-Duty Push Bar	
	DOT TRYKE (12" wheels, fixed of	drive) - 50-FC-1412 Arm L	
	Seating System: Saddle Seat		
•	s: \square Bench Seat \square Pommel S	•	
			push grip, backpad, saddle seat & 2 laterals)
☐ Option #2 — Alternate Seat Bottom Alternates		eries Seat Back System (p Saddle Seat 🗆 Medium (11.5x:	
			x10") □ Skinny Saddle Seat (8.5x6") it □ ProSeries Lumbar Pad
Trunk Support Options	s: ProSeries Full Padded Back	 ☐ Proseries Head Rest ☐ Recumbent 10 Degre 	
	Seating System: Gray Bucke		e Seat Post
•			
•	e Seating System: Snappy Sea	-	
Snappy Seat Accessories:	□ Notched Seat Bottom Cushi		☐ Head Rest
	□ Notched Wide Bottom Cush		Back** ☐ Separator Cube
1412 Accessories:	□ Pull/Steering Bar	☐ Toe Pedal Leveler Pul	•
	☐ Rear Steering Kit		(cannot be used with Rear Steering)
	□ Push Bar		□ 5" (BMX) □ 8" □ Exercise Pedals
	☐ Heavy-Duty Push Bar	☐ Knee Separator: ☐ 2"	□ 5" □ 9" Extender Tube
*cannot be used in conjunction wit	th rear steering. **cannot be used in co	conjunction with an h-harness, head	I rest or laterals.
NFWI High Funct	tioning Set Up (Both Options	c Allow Coasting) Can be us	rad with 1416 1420 1420XI
_	ility to steer, apply brake and mak		
Choose Tryke:		Re a complete pedal rotation ☐ 1420XL	on their own.,
-			· \ A dela the ability to coast and trail ride
		·	ety) Adds the ability to coast and trail ride
		·	ars for hilly terrain — Replaces the 1600 Series
<u> </u>	m (for either option): Saddle Se	·	
Seat Bottom Alternates:	:: □ Bench Seat □ Tractor Seat □ Medium Saddle (10.5x10")		el Saddle Seat □ Med (11.5x10") □ Large(14x13") addle (13x12")
Cuggested Assessaries			

Recipient's Name:			
FOOT TRYKES			
☐ 1416 ProSeries F	OOT TRYKE (16" wheels, fixed	or freewheel drive) – 50-FC-1416 Arr	m Length 14-22", Leg Length 24-32"
□ Option #1 – Standard	Seating System: 1400 ProSer	ries Seat Back System (push grip, b with bracket □ Pommel Saddle So	packpad, saddle seat & 2 laterals)
Trunk Support Options	: □ ProSeries Full Padded Back	□ ProSeries Head Rest□ Recumbent 10 Degree Seat Pos	□ ProSeries Lumbar Pad st
·	e Seating System: Saddle Seat, : □ Bench Seat □ Tractor Seat	/1600 Simple Seat Back with bracket □ Pommel Saddle	Seat □ Med (11.5x10") □ Large(14x13")
☐ Option #3 – Alternate	Seating System: Gray Bucket	Seat	
1416 Accessories:	□ Pull/Steering Bar□ Rear Steering Kit□ Push Bar□ Heavy-Duty Push Bar	 □ Toe Pedal Leveler Pulley □ Calf & Leg Supports □ Hi-Rise Handlebars: □ 8" □ 12" □ Knee Separator: □ 2" □ 5" 	☐ H-Harness: ☐ 8" ☐ 11.5"☐ Exercise Pedals☐ 3-Speed Kit☐ 9" extender tube
1430 Buo Covino E	OOT TRYKE (
		or freewheel drive) – 50-FC-1420 Arr	
☐ Option #1 — Standard Seat Bottom Alternates: Trunk Support Options:	☐ Bench Seat ☐ Tractor S	ies Seat Back System (push grip, ba Seat w/ bracket □ Pommel Saddle □ ProSeries Full Padded Back □ Recumbent 10 Degree Seat Po	Seat □ Med (11.5x10") □ Large(14x13") □ ProSeries Lumbar Pad
☐ Option #2 — Alternate Seat Bottom Alternates	e Seating System: Saddle Seat, : □ Bench Seat □ Tractor S	/1600 Simple Seatback eat w/ bracket □ Pommel Saddle	e Seat □ Med 11.5×10") □ Large(14×13")
☐ Option #3 – Alternate	Seating System: Gray Bucket	Seat	
1420 Accessories:	☐ Pull/Steering Bar	☐ Toe Pedal Pulley	☐ H-Harness: ☐ 8" ☐11.5"
	☐ Rear Steering Kit	☐ Calf & Leg Supports	☐ Exercise Pedals
	☐ Push Bar	☐ 12" Hi-Rise Handlebars	☐ 3-Speed Kit
	☐ Heavy-Duty Push Bar	☐ Knee Separator: ☐ 2" ☐ 5"	□ 9" extender tube
☐ 1420 <u>XL</u> ProSeries	•	xed or freewheel drive) – 50-FC-1420 - 1-32", Leg Length 28-43"	XL
☐ Option #1 – Standard	Seating System: 1400 ProSer	ies Seat Back System (push grip, b	ackpad, saddle seat & 2 laterals)
	ates: □ Bench Seat □ Tractor ions: □ ProSeries Head Rest	Seat w/ bracket □ Pommel Saddle □ ProSeries Full Padded Back □ Recumbent 10 Degree Seat Po	☐ ProSeries Lumbar Pad
☐ Option #2 – Alternate	Seating System: Saddle Seat		
•		or Seat with bracket	☐ Pommel Saddle Seat ¬ ☐ Med (11.5x10") ☐ Large(14x13"
1	e Seating System: Wheelchair sories: Wheelchair Seat Bar Er		
1420 XL Accessories:	☐ Pull/Steering Bar	☐ Toe Pedal Pulley	☐ H-Harness: ☐ 8" ☐11.5"
	□Rear Steering Kit	☐ Calf & Leg Supports	☐ Exercise Pedals
	□ Push Bar	☐ Hi-Rise Handlebars ☐ 12" ☐ 23"	☐ 3-Speed Kit
	☐ Heavy-Duty Push Bar	☐ Knee Separator: ☐ 2″ ☐ 5″	□ 9" extender tube

Recipient's Name:	
OOT TRYKES	
□ NEW! 2722 FOOT TR	RYKE (with 24" wheels, single speed, hand & coaster brake) — 50-FC-2722 Arm Length 22-30", Leg
	Only; must be able to make a pedal revolution on their own
□ Option #1 – Standard Sea	ting System: Large Saddle Seat (13x12")
Seat Bottom Alternate	s: □ Bench Seat □ Med Saddle (10.5x10) □ Tractor Seat w/Bracket: □19″ □23″
	☐ Pommel Saddle Seat ☐ Medium (11.5x10") ☐ Large(14x13")
□ Option #2 – Alternate Sea	ating System: 1400 ProSeries Seat Back (push grip, backpad, saddle seat (13x12") & 2 laterals)
Seat Bottom Option	s: □ Bench Seat □ Med Saddle (10.5x10) □ Tractor Seat w/Bracket: □19" □23"
	□ Pommel Saddle Seat □ Medium (11.5x10") □ Large(14x13")
☐ Option #3 – Alternate Sea	ating System: Large Saddle Seat (13x12")/ 1600 Simple Seat Back
Seat Bottom Option	s: □ Bench Seat □ Med Saddle (10.5x10) □ Tractor Seat w/Bracket: □19" □23"
	☐ Pommel Saddle Seat ☐ Medium (11.5x10") ☐ Large(14x13")
2722 Accessories:	☐ Exercise Pedals ☐ Expanding Pedals ☐ Heel Pedal Leveler Pulley (only w/ Expanding Pedals)
	□ 23" Hi-Rise Handlebars □ Calf & Leg Supports (only w/ Expanding Pedals)
	□ ProSeries Handlebars □ Stem Riser □ 7-Speed Kit □ Dual Hand Brake Lever
☐ JT-2000 FOOT TRYKE	(Recumbent 14-Speed) – 50-FC-2000 (advanced riders only) Arm Length 20-28", Leg Length 30-41"
JT-2000 Accessories:	☐ Exercise Pedals ☐ XL Exercise Pedals ☐ Expanding Pedals ☐ Digital Speedometer
	□ Dual Hand Brake □ Toe Clips □ Basket □ Rearview Mirror
□ JT-2300-USS FOOT TR	YKE (Recumbent 14-Speed with Under Seat Steering) – 50-FC-2300-USS (advanced riders only)
	Arm Length 20-28", Leg Length 30-41"
JT-2300-USS Accessories:	☐ Exercise Pedals ☐ XL Exercise Pedals ☐ Expanding Pedals ☐ Digital Speedometer
	□ Dual Hand Brake □ Toe Clips □ Basket □ Rearview Mirror
	·
Available Summer 2016	;
☐ TADPOLE FOOT TRYK	E (Recumbent 14-Speed) – 50-FC-3000 (advanced riders only) Arm Length 20-28", Leg Length 30-41"
Tadpole Accessories:	□ Exercise Pedals □ XL Exercise Pedals □ Expanding Pedals □ Digital Speedometer
Trapore Accessories.	□ Dual Hand Brake □ Toe Clips □ Rearview Mirror
Therapist Assembly Notes o	r Comments (if any):
The apise hissemal, historica	

HAND-FOOT TRYKE	ES .		
□ New! AM-10 HAN	D-FOOT TRYKE (12" wheels & 3" cran	ık arms) 50-HFC-0105 Arm len	gth 12-16". leg length 14-21"
	Seating System: Blue Bucket Seat	·	, , ,
·	Seating System: Snappy Seat System	 }	
	ories: ☐ Notched Seat Bottom Cushion*		☐ Head Rest
	□ Notched Wide Bottom Cushion'		
AM-10 Accessories:	☐ Rear Steering Kit	☐ Pull/Steering Bar	☐ H-Harness: ☐ 8" ☐ 11.5"
	☐ Push Bar	☐ Heavy-Duty Push Bar	☐ Cube Separator
	D-FOOT TRYKE (12" wheels & 3" crank	arms) 50-HFC-0110 Arm length	12-16", leg length 16-21"
Tryke comes standard w			
•	Seating System: Bench Seat & 1600 S	•	
	Seating System: Snappy Seat System		
Snappy Seat Accesso	ries: ☐ Notched Seat Bottom Cushion*	☐ Laterals	☐ Head Rest
	□ Notched Wide Bottom Cushion*	☐ XL Seat Back**	
	eating System: Blue Bucket Seat		
AM-12S Accessories:	☐ Pull/Steering Bar	☐ Toe Pedal Leveler Pulley	☐ H-Harness: ☐ 8" ☐ 11.5"
		☐ Cube Separator	☐ Exercise pedals
Tryke comes standard w			, leg length 21-29"
Seat Bottom Alternates:	Seating System: Saddle Seat/1600 Sir	•	ledium (11.5x10") □ Large(14x13")
	Seating System: ProSeries Seat Back		
Seat Bottom Alternates:	= :	-	ledium (11.5x10") Large(14x13")
Trunk Support Upgrades:		□ ProSeries Head Rest □	
	Seating System: Gray Bucket Seat		
	Seating System: Snappy Seat System	1	
Snappy Seat Accessories:	☐ Notched Seat Bottom Cushion*		☐ Head Rest
	☐ Notched Wide Bottom Cushion*	☐ XL Seat Back**	☐ Cube Separator
AM-12 Accessories:	☐ Pull/Steering Bar ☐	Toe Pedal Leveler Pulley	☐ H-Harness: ☐ 8" ☐ 11.5"
	☐ Knee Separator: ☐ 2" ☐ 5"		☐ Exercise pedals
	TTRYKE (16" wheels & 5" crank arms) 5	0-HFC-0411 Arm length 18-27	", leg length 24-36"
Now comes standard w	<u> </u>		
☐ Option #1 — Standard : Seat Bottom Alternates:	Seating System: Saddle Seat/1600 Sir		1- di (1(
			ledium (11.5x10") Large(14x13")
	Seating System: 1400 ProSeries Seat		
Seat Bottom Upgrades: Trunk Support Upgrades	☐ Bench Seat ☐ Tractor Seat ☐ ProSeries Full Padded Back	☐ ProSeries Head Rest ☐	ledium (11.5x10") □ Large(14x13") ProSeries Lumbar Pad
	Seating System: Gray Bucket Seat	_ 1103CHC3 HCdd Ne3t _	1 100cmed Eumbur 1 au
AM-16 Accessories:		edal Leveler Pulley	☐ H-Harness: ☐ 8" ☐ 11.5"
10 / 10003301103.		•	e pedals

Recipient's Name:

^{*}cannot be used in conjunction with rear steering. **cannot be used in conjunction with an h-harness, head rest or laterals.

ke is currently being u	used for research!		0 Arm length 18-30", leg length 32-46"
1-20 Accessories:	☐ Exercise Pedals	☐ XL Exercise Pedals	☐ Expanding Pedals
	☐ Wheelchair Seat Bar Ends	☐ Swing Away Arms	☐ Calf & Leg Supports (Expanding ped
eranist Assembly Note	s or Comments (if any):		
erapiser issembly riote	o or comments (ii any).		

□ New! AM-10 HAND	TRYKE (12" wheels, 2.5" crank arms &	permanent foot platform)	50-HFC-0105 Arm length 12-16"
	eating System: Blue Bucket Seat		
□ Option #2 – Upgraded S	eating System: Snappy Seat System		
Snappy Seat Accesso	ries: ☐ Notched Seat Bottom Cushion*	□ Laterals	☐ Head Rest
	☐ Notched Wide Bottom Cushion	* □ XL Seat Back**	
AM-10 Accessories:	□ Rear Steering Kit	☐ Pull/Steering Bar	☐ H-Harness: ☐ 8" ☐ 11.5"
	□ Push Bar	☐ Heavy-Duty Push B	ar □ Separator Cube
☐ AM-12Small HAND	TRYKE (12" wheels, 3" crank arms & perr	manent foot platform) 50 -	HFC-0110 Arm length 12-16"
Tryke comes standard wit	th rear steering!		
□ Option #1 – Standard Se	eating System: Bench Seat & 1600 Si	mple Seat Back	
	eating System: Snappy Seat System		
Snappy Seat Accesso	ries: ☐ Notched Seat Bottom Cushion*		☐ Head Rest
	□ Notched Wide Bottom Cushion	* □ XL Seat Back**	
•	eating System: Blue Bucket Seat	-	
AM-12S Accessories:	☐ Pull/Steering Bar ☐	•	☐ H-Harness: ☐ 8″ ☐ 11.5″
		L	Separator Cube
	– (12" wheels, 4" crank arms & permane	nt foot platform) 50-HC-02	210 Arm Length 14-23"
Tryke comes standard wit		vala Caat Daal	
Seat Bottom Alternat	eating System: Saddle Seat/1600 Sim		t □ Medium (11.5x10") □ Large(14x13")
	eating System: 1400 ProSeries Seat I		
	nates: Bench Seat		t □ Medium (11.5x10") □ Large(14x13")
	tions: ProSeries Full Padded Back		t □ ProSeries Lumbar Pad
	eating System: Gray Bucket Seat		
•	eating System: Snappy Seat System		
	ries: Notched Seat Bottom Cushion*	□ Laterals	☐ Head Rest
	☐ Notched Wide Bottom Cushion	* □ XL Seat Back**	☐ Separator Cube
AM-12 Accessories:	☐ Pull/Steering Bar	☐ Toe Pedal Leveler	Pulley ☐ H-Harness: ☐ 8" ☐ 11.5"
	☐ Knee Separator: ☐ 2" ☐ 5"	☐ 9" Extender Tube	
cannot be used in conjunction with rea	r steering. **cannot be used in conjunction with	n an h-harness, head rest or later	als.
☐ AM-16 HAND TRYKE	(16" wheels, 5" crank arms & permanent t	oot platform) 50-HFC-04:	11 Arm length 18-27"
Tryke comes standard wit	h rear steering!		
•	ating System: Saddle Seat/1600 Sim		
Seat Bottom Alternates:			t □ Medium (11.5x10") □ Large(14x13")
	eating System: 1400 ProSeries Seat	•	
	☐ Bench Seat ☐ Padded Tractor Seat		_
	☐ ProSeries Full Padded Back	☐ Proseries Head Resi	t □ ProSeries Lumbar Pad
¬ option #3 − opgraded S	eating System: Gray Bucket Seat		
AM-16 Accessories:	☐ Pull/Steering Bar ☐	Toe Pedal Leveler Pulley	☐ H-Harness: ☐ 8" ☐ 11.5"

Recipient's Name:

Recipient's Name:		
HAND TRYKES		
□ 1024 HAND TRYKE	fka 1424 – (24" wheels & 3-Spe	ed drive train) 50-HC-1024 Arm length 22-26"
1024 Accessories:	☐ Wheelchair Seat Bar Ends	☐ Swing Away Arms ☐ Remote location hand brake & shifting kit
		ont wheel & 7-Speed grip shift) 50-HC-1000 Arm length 22-26"
HP-1000 Accessories:	☐ Wheelchair Seat Bar Ends	☐ Swing Away Arms
-	•	
Therapist Assembly Note	es or Comments (if any):	

Amtryke Sizing Chart

TRYKE TYPE (How will the tryke be propelled?)	RIDER LEG LENGTH (Inches from center of hip to bottom of shoe.)	RIDER ARM LENGTH (Inches from middle of shoulder to center of digit crease.)	MODEL	RIDER WEIGHT (Pounds)	RIDER MAX HEIGHT (Inches)	TRYKE WEIGHT (Pounds)	WHEEL SIZE (Inches)	TTRYKE HEIGHT (inches)	TRYKE LEGNTH (inches)	TRYKE WIDTH (Inches)		
Hand & Foot	14-21	12-16	AM-10	55	40	45	10	24	38	21		
	19-24	15-20	AM-12S	150	40	45	12	27	38	24		
	21-29	14-23	AM-12	150	47	45	12	36	60	32		
	24-36	18-27	AM-16	175	66	55	16	36	68	33		
	32-46	18-30	AM-20	250	74	80	20	45	72	32		
Foot	9-16	11-20	1410	55	40	45	10	24	38	21		
	17-22	9-13	Snappy	150	42	45	12	30	38	24		
	21-26	14-22	1412	125	42	72	12	13	43	27		
	24-30	14-22	1416	175	60	74	16	49	58	30		
	27-33	18-28	1420	250	68	74	20	50	64	30		
	28-43	20-32	1420XL	275	76	89	20	43	72	29		
	28-36	22-30	2722	275	72	76	24	42	60	30		
	30-41	20-28	JT-2000/ JT2300USS	250	74	80	20	48	72	32		
Hand	22-41	22-26	1024	250	72	74	20	45	75	32		
	22-41	22-26	HP-1000	250	74	80	Front 16 Rear 20	45	72	32		
I	All trykes in the Hand & Foot section can be converted to Hand trykes – except the AM-20.											

A Center of Shoulder


B Center of Elbow

C Center of Digit Crease

D Center of Hip (greater trochanter)

E Center of Knee

F Bottom of Foot



RIDER'S MEASUREMENTS

Arm Measurements (inches) Total Length

Left A to B + B to C = _____

Right A to B + B to C = ______ Trunk A to D = _____

Leg Measurements (inches) Total Length

Left D to E + E to F = _____

Right D to E + E to F = _____

HELMET SIZING

Youth Sizes Head Circumference Inches

XXS 18.5 to 19.5 XS 20.5 to 22

S/M 22 to 23.6 L/XL 23.6 to 25.75

Adult Sizes Head Circumference Inches

S/M 22 to 23.6 L/XL 23.6 to 25.75